


## ESTi™ Analysis - Molly

Muscles

Ligaments


| 0 sec | | 5 sec | | 10 sec | | 15 sec | | 20 sec | | 25 sec | | 30 sec | | 35 sec | | 40 sec | | 45 sec | | 5

Left 1



E	S	T	Threshold	Max T	Max S
0.5	6.5	6.2	0.1	250	0.99

Right 1



2.2	7.6	6.5	0.1	250	0.99
-----	-----	-----	-----	-----	------




FIRST: Select your reading frame – avoid the start and finish sequences unless you are interested in them ..  
Just place the cursor on the edge of the green measuring field to expand or shrink it

# ESTi™ Analysis - Molly

Muscles

Ligaments

| 0 sec | | 5 sec | | 10 sec | | 15 sec | | 20 sec | | 25 sec | | 30 sec | | 35 sec | | 40 sec | | 45 sec | | 5



E	S	T	Threshold	Max T	Max S
0.5	6.5	6.2	0.1	250	0.99

2.2 7.6 6.5 0.1 250 0.99

NEXT: Adjust your threshold .. It should be above the background noise and nicely in the data

adjust threshold here

# ESTi™ Analysis - Molly

Muscles

Ligaments

| 0 sec | | 5 sec | | 10 sec | | 15 sec | | 20 sec | | 25 sec | | 30 sec | | 35 sec | | 40 sec | | 45 sec | | 5

Left 1

Right 1

E	S	T	Threshold	Max T	Max S
0.5	6.5	6.2	0.1	250	0.99

2.2	7.6	6.5	0.1	250	0.99
-----	-----	-----	-----	-----	------

Adjust T-max if needed (250 = 250Hz) – if you have a fast contracting muscle then the pre-set T-max may be too low and will therefore give you 0 values if not raised


# ESTi™ Analysis - Molly

Muscles


Ligaments

| 0 sec | | 5 sec | | 10 sec | | 15 sec | | 20 sec | | 25 sec | | 30 sec | | 35 sec | | 40 sec | | 45 sec | | 5

Left 1



Right 1



E	S	T	Threshold	Max T	Max S
0.5	6.5	6.2	0.1	250	0.99

2.2	7.6	6.5	0.1	250	0.99
-----	-----	-----	-----	-----	------




You can adjust the S-max too if you have a small amplitude signal .. this just alters the amplitude scale 0.99 = 1 volt, so 0.50 = 0.50 volt etc


# ESTi™ Analysis - Molly

Muscles

Ligaments



Now you can read your E, S & T-scores for the left side and right side. The S and T scores are inverted .. that is to say that a large amplitude signal will give a low S-score and a high firing frequency a low T-score



For a signal duration of  $105(E1)+120(E2)+75(E3) = 300$  mSec over a recording timeframe (TT) of 900 mSec the E-score would be:

$$E = \frac{900-300}{900} * 10 = 6.6$$

For an amplitude of 450 mV in relation to a full 999 mV signal (6 dB), the S-score would be:

$$S = \frac{999-450}{999} * 10 = 5.5$$

For a signal frequency of 53 Hz (16 spikes per 300 mSec) in relation to a maximum level of 120 Hz, the T-score would be:

$$T = \frac{120-53}{120} * 10 = 5.6$$

Thus the overall ESTi-score would be calculated as:

$$ESTi = \frac{6.6+5.5+5.6}{3} * 10 = 5.9$$


Here is how the E, S and T scores are calculated using the amplitude and the frequency settings (S-max and T-max)

The E-score is calculated knowing the total recording frame time and calculating the periods of time the muscle is active during that total time.

# ESTi™ Analysis - Molly

Muscles

Ligaments



Using the Left E,S,T and the Right E,S,T scores you can calculate the muscle balance ...

$((0.5+6.5+6.2)-(2.2-7.6-6.5)) = 13.2-16.3 = -3.1$  (this means the left side is working slightly more than the right side)

# ESTi™ Analysis - Molly

Muscles

Ligaments

| 0 sec | | 5 sec | | 10 sec | | 15 sec | | 20 sec | | 25 sec | | 30 sec | | 35 sec | | 40 sec | | 45 sec | | 5

Left 1



Right 1

E	S	T	Threshold	Max T	Max S
0.5	6.5	6.2	0.1	250	0.99

2.2	7.6	6.5	0.1	250	0.99
-----	-----	-----	-----	-----	------



Hope this helps .. have fun 😊