

PRODUCT DATA SHEET

TECHNICAL SPECIFICATIONS

Name of the unit:	CURO
Number of channels:	4
Memory:	4GB (SD card)
Built-in memory back-up:	EEPROM
Sampling rate:	2 kbps
Operating time:	4 hours
Power supply:	12V/1A
Dimensions:	W: 9.3, L: 14.4, H: 3.5 cm
Weight:	318 g

ACCESSORIES

AMG preamplifier for CURO

Number of gain settings:	2 (6dB/15dB)
Sampling rate:	2 kbps
Power supply:	Powered by CURO
Length of wire:	2.1 meter
Type of connector:	RJ45

AMG surface electrodes

Type:	Piezoelectric crystals
Sizes:	Ø 2 cm and Ø 5 cm
Length of wire:	8.2 cm
Type of connector:	Jack

Data transfer cables

Type:	mini USB to USB
Length:	2.0 meter

Recharger cable

Type:	AC-DC,
Output:	12V/0.66A
Length:	1.8 meter

